


# CLEAN YOUR HANDS, Per Facility Guidelines, During The Disinfection Cycle.

## Easy to Use

- 1 Clean visible dirt, moisture or liquid residue from items to be disinfected as per manufacturer's instructions



- 2 Make sure the unit is on -the green light will be lit



- 3 Place items inside without stacking



- 4 Close the door and turn the lock knob



- 5 Press the red button to begin disinfection

**Caution: Place and Remove Items Carefully.**

**AUVS**  
ADVANCED ULTRA-VIOLET SYSTEMS

For Sales and Service: Safety Net, LLC  
866-918-8546 | [info@safetynetamerica.com](mailto:info@safetynetamerica.com)

**Safety Net**  
Solutions that save lives.

